CHOLESTEROL HDL/LDL CALIBRATOR

REF 5030301 3x1 ml
IVD For in vitro medical device
Per uso diagnostico in vitro

Uso

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Use

This product is intended for in vitro use in the calibration of HDL and LDL cholesterol clearance methods.

Summary

Lyophilized calibrator serum prepared from human serum.
The concentration and activities of the calibrator components have been adjusted for a wide range of automatic analyzers to ensure optimal calibration of the clinical chemistry reagents.

Calibration values

The concentration and activities of the components are lot-specific.

LOT N° 1728CH

ANALYTE	SYSTEM	TARGET	
		mmol/l	mg/dl
Direct HDL	Mean of all instruments	1.91	73.7
Direct LDL		4.33	167

Preparation

Open one bottle of calibrator avoiding the loss of lyophilizate.
Reconstitute the lyophilized serum with exactly 1.0 ml of bi-distilled water.
Close the bottle and let it stand for 30 minutes out of bright light.
Before to use mix by inverting the vial.
Do not shake the vial as the formation of foam should avoid.

Storage and stability

Prior to reconstitution, when stored at 2-8 °C, the calibrator is stable until the expire date.

After the reconstitution the serum is stable:

5 days	if stored at	+4° C
30 days	if stored at	-20° C

Note

- Avoid direct light exposure, contamination and evaporation.
- In case of compliant or quality control request, refer to the lot number on the package or to the lot number of the single vial.
- Non properly recostituion and / or storage can feed erratic results.

Precaution in use

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998).

Human serum is obtained using only blood of donors tested by FDA method and found non-reactive for HbsAg and negative for antibodies to HIV-1/2 and HCV. However, no known test method can offer complete assurance that product derived from Human blood transmit infectious diseases, this product should be handled as potentially infectious biological material with caution according to good laboratory practice.

Waste management

Please refer to the local legal requirements.

Bibliography

Department of Labor, Occupational Safety and Health Administration 29 CFR Part 1910.1030 Fed. Register July 1, 1998; 6:267-280
Council Directive 90/679/EEC. Official journal of the European directives N°L374 from Dec.31,1990:1-12

Symbols

	Biological risk
	CE Mark (98/79 CE regulation)
	in vitro medical device
	Batch Code
	Use by
	Storage temperature limits
	Read instruction for use
	Gesan Production srl